

Painted Steel Cast Cart/ Model: 6030

Eight Drawer, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Stainless steel bucket holder with vinyl mat
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free
- Designed in conjunction with orthopedic technicians to maximize functionality

Other Cast Cart Options:

- All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- 43" H x 38" W x 21" D (109.2 cm H x 96.5 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 171 lbs (77.7 kg)
- Approximate shipping weight: 216 lbs (98.2 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)