

## Painted Steel Cast Cart/ Model: 6035

Eight Drawer, Deluxe Package



### Seventeen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

### Construction:

- Painted steel cabinet and aluminum drawers

### Features:

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet/ pump pole
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free

### Other Cast Cart Options:

- All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

### Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

### Overall Cart Dimensions:

- 43" H x 38" W x 21" D  
(109.2 cm H x 96.5 cm W x 53.3 cm D)  
w/ pump pole 76" H x 44" W x 21" D  
(193 cm H x 111.8 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 202 lbs (91.8 kg)
- Approximate shipping weight: 247 lbs (112.3 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D  
(144.8 cm H x 119.4 cm W x 68.6 cm D)