

EL51LT



33.5" x 67" x 27.25" (H x W x D)

Commercial -45° C capable chest freezer with digital thermostat and 15.5 cu.ft. capacity

Highlights:

Capable of -45° C operation for select laboratory and industrial storage purposes

UL-S listed to NSF-7 commercial standards

Factory installed lock for security you can depend on

Product Features:

Commercially approved	UL-S listed to NSF/ANSI-7 standards
-45° capable	Low temperature operation is ideal for a variety of laboratory, life science, and industrial applications
Extra thick insulation	4 1/2" thick walls for added efficiency and stable cooling
Factory installed lock	Keyed lock on the handle for reliable security
Pre-installed probe hole	Factory-installed 1/2" probe hole located on the right side to accept user-provided monitoring equipment
Digital thermostat	Electronic controls for precise temperature management, with temperature readout and controls externally located for added convenience
Counterbalanced hinges	Easy and safe loading and unloading with a lid that remains open when lifted past 45° angles
Casters	Four solid wheels included for optional mobility
Large capacity	Over 15 cu.ft. of interior storage space

EL51LT Specifications:

Overview	
Height of Cabinet	33.5" (85 cm)
Width	67.0" (170 cm)
Depth	27.25" (69 cm)
Depth with Handle	28.75" (73 cm)
Capacity	15.5 cu.ft. (439 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	7.1
Voltage/Frequency	115 V AC/60 Hz
Weight	210.0 lbs. (95 kg)
Shipping Weight	220.0 lbs. (100 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101023519
Freezer	
Thermostat Type	Digital
Refrigerant Type	R290a
Refrigerant Amount	3.35 oz.
High Side PSI	342.0
Low Side PSI	174.0
Wheels	4

