

FF6LWPLUS2ADA



32.25" x 23.63" x 23.5" (H x W x D)

24" wide freestanding ADA height all-refrigerator in white, auto defrost with a lock, NIST calibrated thermometer, digital thermostat, door storage, and internal fan

Highlights:

32" height for use in ADA compliant settings

Ideal 24" footprint offers large capacity in a slim fit

Factory installed lock provides security you can count on

Product Features:

| | |
|-------------------------------------|---|
| ADA compliant | Sized at 32" high for use in ADA compliant settings |
| Automatic defrost | Reduced maintenance with auto defrost system |
| NIST calibrated temperature display | Thermometer provides continuous display of the current and high/low temperature to the nearest tenth of a degree, with a readout switchable from Celsius to Fahrenheit at the touch of a button |
| Buffered temperature probe | Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature |
| Digital thermostat | Electronic controls ensure easy and intuitive temperature management |
| Factory installed lock | Keyed lock for secure interior |
| Internal fan | Fan-forced cooling ensures the unit is evenly cooled from top to bottom |
| Adjustable shelves | Rearrange refrigerator space to accommodate all sizes or remove shelves for simple cleanup |
| Door storage | Store items right on the door for easy access |

FF6LWPLUS2ADA Specifications:

| Overview | |
|----------------------------|-----------------------|
| Height of Cabinet | 32.0" (81 cm) |
| Height to Hinge Cap | 32.25" (82 cm) |
| Width | 23.63" (60 cm) |
| Depth | 23.5" (60 cm) |
| Depth with Handle | 25.13" (64 cm) |
| Depth with door at 90° | 46.13" (117 cm) |
| Capacity | 5.5 cu.ft. (156 L) |
| Defrost Type | Automatic |
| Door | White |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 1.4 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 100.0 lbs. (45 kg) |
| Shipping Weight | 105.0 lbs. (48 kg) |
| Parts & Labor Warranty | 1 Year |
| Compressor Warranty | 5 Years |
| Refrigerator Features | |
| Door Swing | RHD |
| Reversible | Factory Reversible |
| Shelf Type | Wire |
| Shelf Qty | 3 |
| Full Door Shelf Qty | 3 |
| Adjustable Shelves | Yes |
| Thermostat Type | Digital |
| Fan Type | Interior |
| Refrigerant Type | R600a |
| Refrigerant Amount | 0.64oz. |
| High Side PSI | 180.0 |
| Low Side PSI | 80.0 |
| Level Legs Qty | 4 |
| Interior Light | Yes |
| Temperature Range | 2 to 8°C / 36 to 46°F |
| Dimensions | |
| Interior Height | 27.75" (70 cm) |
| Interior Width | 20.0" (51 cm) |

| | |
|------------------------|---------------|
| Interior Depth | 17.5" (44 cm) |
| Compressor Step Height | 7.0" (18 cm) |
| Compressor Step Width | 20.0" (51 cm) |
| Compressor Step Depth | 6.0" (15 cm) |

