

# FF7LWPLUS2



33.5" x 23.63" x 23.5" (H x W x D)

24" wide freestanding all-refrigerator in white, auto defrost with a lock, NIST calibrated thermometer, digital thermostat, and internal fan

## Highlights:

Ideal 24" footprint offers large capacity in a slim fit

Fan-forced cooling ensures the unit is evenly cooled from top to bottom

Thermometer provides external readout of current and high/low temperature

## Product Features:

|                                     |                                                                                                                                                                                                 |
|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Commercially approved               | ETL-S listed to NSF standards for commercial use                                                                                                                                                |
| NIST calibrated temperature display | Thermometer provides continuous display of the current and high/low temperature to the nearest tenth of a degree, with a readout switchable from Celsius to Fahrenheit at the touch of a button |
| Buffered temperature probe          | Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature                                            |
| Digital thermostat                  | Electronic controls ensure easy and intuitive temperature management                                                                                                                            |
| Factory installed lock              | Keyed security on the refrigerator door                                                                                                                                                         |
| Internal fan                        | Fan-forced cooling ensures the unit is evenly cooled from top to bottom                                                                                                                         |
| Automatic defrost                   | Reduced maintenance with auto defrost system                                                                                                                                                    |
| Adjustable shelves                  | Rearrange refrigerator space to accommodate all sizes or remove shelves for simple cleanup                                                                                                      |
| Flat door liner                     | Easy to clean, with more depth for storage                                                                                                                                                      |

## FF7LWPLUS2 Specifications:

| Overview                   |                       |
|----------------------------|-----------------------|
| Height of Cabinet          | 33.5" (85 cm)         |
| Width                      | 23.63" (60 cm)        |
| Depth                      | 23.5" (60 cm)         |
| Depth with Handle          | 25.13" (64 cm)        |
| Depth with door at 90°     | 46.13" (117 cm)       |
| Capacity                   | 5.5 cu.ft. (156 L)    |
| Defrost Type               | Automatic             |
| Door                       | White                 |
| Cabinet                    | White                 |
| US Electrical Safety       | ETL                   |
| Canadian Electrical Safety | ETL-C                 |
| Sanitation                 | ETL-S                 |
| Amps                       | 1.4                   |
| Voltage/Frequency          | 115 V AC/60 Hz        |
| Weight                     | 100.0 lbs. (45 kg)    |
| Shipping Weight            | 105.0 lbs. (48 kg)    |
| Parts & Labor Warranty     | 1 Year                |
| Compressor Warranty        | 5 Years               |
| Refrigerator Features      |                       |
| Door Swing                 | RHD                   |
| Reversible                 | Factory Reversible    |
| Shelf Type                 | Wire                  |
| Shelf Qty                  | 3                     |
| Adjustable Shelves         | Yes                   |
| Thermostat Type            | Digital               |
| Fan Type                   | Interior              |
| Refrigerant Type           | R600a                 |
| Refrigerant Amount         | 0.64oz.               |
| High Side PSI              | 180.0                 |
| Low Side PSI               | 80.0                  |
| Level Legs Qty             | 2                     |
| Interior Light             | Yes                   |
| Temperature Range          | 2 to 8°C / 36 to 46°F |
| Dimensions                 |                       |
| Interior Height            | 27.75" (70 cm)        |
| Interior Width             | 20.0" (51 cm)         |
| Interior Depth             | 17.5" (44 cm)         |

|                        |               |
|------------------------|---------------|
| Compressor Step Height | 7.0" (18 cm)  |
| Compressor Step Width  | 20.0" (51 cm) |
| Compressor Step Depth  | 6.0" (15 cm)  |

