

VLT650



33.5" x 23.63" x 23.5" (H x W x D)

Built-in undercounter laboratory freezer capable of -35° C (-31° F) operation

Highlights:

Commercially listed to NSF-7 standards

Capable of operating at -35° C for long-term food preservation and safe storage of pharmaceuticals

Slim width and counter height fits easily where space is tight

Product Features:

Commercially approved	Meets NSF-7 commercial standards
Slim counter height dimensions	Less than 24 inches wide to fit tight spaces
-35° C (-31° F) capable	Low temperature operation ideal for vaccines, pharmaceuticals, and long-term preservation of frozen food
Built-in capable	Front-breathing system allows the unit to be used built-in
Factory installed lock	Keyed lock offers convenient security
Three slide-out drawers	Conveniently organize your items while keeping the temperature stable when the door is opened
Digital thermostat	Electronic controls for precise temperature management
Alarm with temperature display	Audible alarm sounds to warn of temperature variations with a digital readout on the thermometer to provide the interior temperature
Manual defrost	Static cooling system with manual defrost operation

VLT650 Specifications:

Overview	
Height of Cabinet	33.5" (85 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Capacity	3.5 cu.ft. (99 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Voltage/Frequency	115 V AC/60 Hz
Shipping Weight	110.0 lbs. (50 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101009278
Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	2.76 oz.
Interior Drawers	Yes
Drawers Qty	3

