



# M-series

## MDS1824E05

### M-SERIES SHORT MINI WIDTH CART

Five Drawers with an Electronic Pushbutton Lock

#### M-Series Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

#### Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

#### Colors:

Shown in White with Navy drawers, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).



Engineered & Manufactured in the U.S.



Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 36.125" (91.8 cm)
- Cart dims: 36.38" H x 23.13" W x 22" D (92.4 cm H x 58.7 cm W x 55.9cm D)
- Cart net weight: 86 lbs (39 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300