



# M-series

# MDS1830B06+MD18-EMG

## M-SERIES TALL EMERGENCY CART with MD18-EMG Emergency Package

### Tall Mini Width, Six Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

### MD18-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



Aluminum 16" Push Side Rail (SIDERAIL16)



IV pole (MD-IV2)



Engineered & Manufactured in the U.S.

### Features:

- Sturdy, uni-body steel construction
- Durable scratch resistant heat sealed powder coat finish
- Breakaway lock bar or your choice of other lock options
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
  - 1 steering caster and 2 braking casters
- Body has pre-threaded accessory mounting holes

### Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Color shown in Red. More than 20 additional colors are available.



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
medicalsales@harloff.com • www.harloff.com

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)
- 12.5" W x 17" D (31.8 cm W x 43.2 cm D)

#### Assembled Cart Dimensions (without accessories):

- 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)