



M-series

MDS1830B06

M-SERIES TALL MINI WIDTH CART with Six Drawers and a Breakaway Lock

M-Series Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
 - 1 steering caster and 2 braking casters
- Latex free

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

Colors:

Shown in Red, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.



Engineered &
Manufactured
in the U.S.



Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 41.5" (105.4 cm)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300