



M-series MDS3030B05-12FD

M-SERIES CART WITH FLIP-UP DOOR Crash Cart w/Five Drawers and Breakaway Lock

Includes a lower flip-up door that stows horizontally and, when opened, provides a 12" high inside storage compartment

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- 30" of vertical drawer space
- Cart is delivered assembled
- Latex free



**Color shown is Red.
More than 20 additional
colors are available.**

Tall, Standard Width, Five Drawer Emergency Cart with Flip-Up Door and Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" flip-up drawer (30.4 cm)

Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 150 lbs (68 kg)
- Approximate shipping weight: 189 lbs (85.7 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered &
Manufactured
in the U.S.**



View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300