



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**MRI-Compatible, Six Drawer Anesthesia Cart, Key Lock,
 Specialty Package
 #MR6K-MAN**



Specifications

Construction:

- Non-ferrous painted aluminum cabinet and drawers

Features:

- Designed for compatibility in MRI suites
- Tested for safe use in MRI suites by independent testing authority (copies of report available)
- Cabinet painted gloss white with MRI identification sticker applied
- Over 20 different drawer colors available
- MRI-compatible accessories available
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Corner bumpers, protecting the cart
- Painted aluminum pull-out shelf
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- Twelve year warranty
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (2) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Over 20 drawer colors available:

***Antimicrobial Finish Available**

Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

**Anesthesia Specialty Package Includes
 MRI-Compatible Accessories:**

- Aluminum 3-tier raised back rail system (MR-RBR3)
- Aluminum side rails (4) (MR-SR)
- Plastic full drawer tray organizers (2) (68530-P1)

Assembled Cart Dimensions

- 68.5" H x 32" W x 25" D
 (174 cm H x 81.3 cm W x 63.5 cm D)
- Cart net weight: 137 pounds (62.3 kg)
- Approx. shipping weight: 157 pounds (71.4 kg)