



MRI-SAFE

MR7K-MAN

MR-CONDITIONAL ANESTHESIA CART Seven Drawers, Key Lock & Specialty Package

Specialty Package Includes:



3-Tier Back Rail
(MR-RBR3)



(4) 12" Side Rails
(MR-SR)



(2) Plastic Trays w/
Adjustable Dividers
(68530-P1)



Engineered &
Manufactured
in the U.S.

Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Protective corner bumpers
- Pull-out shelf
- Antimicrobial plastic top with integrated push handles, removable
- Latex free cart is delivered fully assembled with pre-threaded accessory mounting holes
- Durable, heat sealed powder coat finish

Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)
- One 4" drawer (10.2 cm), four 3.25" drawers (8.3 cm), one 6.5" drawer (16.5 cm) & one 9.75" drawer (24.8 cm)



Color shown
with a White body
and Navy drawers.
More than 20
additional colors
are available.

Assembled Cart Dimensions:

- Writing surface: 43.5" (110.5cm)
- 46" H x 31.8" W x 23.75" D (116.8 cm H x 80.8 cm W x 60.3 cm D)
- Cart net weight: 127 lbs (57.6 kg)
- Approximate shipping weight: 166 lbs (75.3 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Specially Made for Magnetic Resonance Imaging Suites

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300