

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

MR-Conditional Line/ Model: MR7K-MAN

Seven Drawer Anesthesia Cart, Key Lock, Specialty Package



Over 20 Drawer Colors Available:

Antimicrobial Finish Available: Lemon Drop Yellow Only

BBurgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

 Non-magnetic or "weekly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed for compatibility in MRI suites
- Tested by independent testing authority to Tesla 3 level (copies of report available)
- Meets American Society for testing and materials (ASTM) international, designation: F2503-5 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel one with brake, one directional
- BEST® brand lock on cabinet, with two keys
- Corner bumpers, protecting the cart
- Painted aluminum pull-out shelf
- Antimicrobial plastic top, removable
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat finish

Anesthesia Specialty Package Includes MRI-Compatible Accessories:

- Aluminum 3-tier raised back rail system (MR-RBR3)
- Aluminum side rails (4 ea.) (MR-SR)
- Plastic full drawer tray organizer (2 ea.) (68530-P1)

Drawer Configuration:

- (1) 4" (10.2 cm)
- (4) 3" (8.3 cm)
- (1) 6" (16.5 cm)
- (1) 9" (24.8 cm)

Internal Drawer Configuration:

• 23.9" W x 17" D (60.8 cm x 43.2 cm)

Overall Cart Dimensions:

- 68.2" H x 32" W x 24" D (173.4 cm H x 81.3 cm W x 61 cm D)
- Cart net weight: 127 pounds (57.6 kg)
- Approximate shipping weight: 166 pounds (75.3 kg)
 Approximate shipping dims: 55" H x 37" W x 26" D
- Approximate snipping dims: 55 H x 37 W x 26 L (139.7 cm H x 94 cm W x 66 cm D)