



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

V-Series, Four Drawer Crash Cart, Specialty Accessory Package #V24-4EMG

This crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price.

This emergency cart includes all commonly required accessories and is perfect for hospitals, clinics, dialysis centers, doctors' offices, urgent care facilities and nursing homes needing to buy carts to store, transport and dispense critical supplies used for emergency situations.



Specifications

Construction:

- Steel construction

Features:

- Break-away lock bar
- 100 break-away seals included
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year warranty
- Ships with wheels installed and ready to use
- Minimal assembly required at delivery
- Latex free
- Durable, hard-baked powder coat finish

EMERGENCY Package Accessories Include:

- Cardiac board and brackets
- Side mounted drop shelf
- IV pole
- Oxygen tank holder (D/E size cylinders)
- Powder coated steel top rail / push handle

Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Integrated Pull Handle on all drawers

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 36.25" H x 33.25" W x 22" D (92.1 cm H x 84.5 W x 55.9 D) (height includes IV Pole)
- Cart weight: 110 lbs (49.9 kg)
- Approx shipping weight: 132 pounds (59.9 kg)