



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

MR-Conditional, Seven Drawer Anesthesia Cart, Key Lock, Standard Package #MR7K

Specifications



Drawer Configuration:

- (1) - 4" (10.2 cm)
- (4) - 3.25" (8.3 cm)
- (1) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Over Twenty available color choices:

*Antimicrobial Finish Available

Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Construction:

- Non-magnetic or "weakly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested for safe use in MRI suites by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Over 20 different drawer colors available
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – two with brake
- BEST® brand lock on cabinet, with 2 keys
- Corner bumpers protecting the cart
- Painted aluminum pull-out shelf
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish

Internal Drawer Dimensions:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Assembled Cart Dimensions

- 45.75" H x 32" W x 25" D
(116.2 cm H x 81.3 cm W x 63.5 cm D)
- Cart net weight: 137 pounds (62.3 kg)
- Approx. shipping weight: 157 pounds (71.4 kg)