



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## Four Drawer Mini Line Phlebotomy Cart Package with Key Lock #3254KPB



### Features

- Wrap-around, stabilizing bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key-locking system with 2 keys
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D (148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black