

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

OptimAL, Aluminum 294 Capacity Unit Dose Medication Cart, Standard Package #AL294BOX

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray,
Mauve, Red, Teal, Yellow, Forest Green, Hammertone
Blue, Hammertone Gray, Hammertone Green, Hammer
Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

· Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 294 unit dose medication box capacity, single width narcotics drawer
- · BEST® brand lock on cabinet with four keys
- · Illinois lock on narcotics drawer with four keys
- · Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel,
 2 with toe locking brake and 1 directional
- · Hard-baked powder coat finish
- Dividers are available to accommodate different medication packaging

Left Column Drawer Configuration:

- 1 ea 4" multipurpose storage drawer
- 6 ea 3.25" drawers with internal dividers
- 1 ea 9.75" multipurpose storage drawer with adjustable dividers

Right Column Drawer Configuration:

- 1 ea 4" multipurpose storage drawer
- 1 ea 9.75" locking narcotics drawer
- 2 ea 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 162 pounds (73.6 kg)
- Approx. shipping weight: 205 pounds (93.2 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 68.6 cm D)