

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

OptimAL Line, Aluminum 460 Punch Card Medication Cart #AL460PC



Specifications

Construction:

· Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 460 punch card capacity, single width narcotics drawer
- BEST® brand lock on cabinet and narcotics drawer
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Painted steel pull-out shelf
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- · Hard-baked powder coat finish
- Offered in specialty package (#AL460PCSP) which includes #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, AL680436 aluminum waste container with lid and #684810 locking sharps container

Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY) Hammertone Green (-HTGN) and Sand (-SD).

Assembled Cart Dimensions

- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 137 lbs (62.2 kg)
- Approx shipping weight: 197 lbs (89.6 kg)
- Writing surface: 44"

Left Column Drawers:

- 1 ea 4" multipurpose storage drawer
- 3 ea 9.75" punch card drawers with dividers

Right Column Drawers:

- 1 ea 4" multipurpose storage drawer
- 1 ea 9.75" locking narcotics drawer
- 2 ea 9.75" multipurpose storage drawer

Internal Drawer Dimensions:

- Left Side: 17" D x 20.5" W (43.2 cm x 52.07 cm)
- Right Side: 17" D x 6.25" W (43.2 cm x 15.9 cm)