



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Triple Column Mobile Medical Storage Cabinet with Key Locking Roll-up Tambour Door #MS8160-CARDIO

The MS8160-CARDIO is a triple column medical storage cabinet secured by a key locking roll-up tambour door. Utilizing the most innovative modular tray system on the market, this mobile cabinet can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.

The MS8160-CARDIO is preconfigured as follows: left column, top to bottom, 3 ea MS-SH shelves, 2 ea 2" trays, 2 ea 4" trays, 3 ea 8" trays; middle column, top down, 4 ea 2" trays, 3 ea 4" trays, 2 ea 5" wired baskets, 2 ea 8" wired baskets; right column, top to bottom, 3 ea MS-SH shelves, 2 ea 2" trays, 2 ea 4" trays, 3 ea 8" trays



### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- Maximize space and minimize costs using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage, reducing capital purchase cost.
- Pull-out-stoppers available for trays, allowing trays to be viewed at an angle without falling out of the cart.
- Catheter shelf can be added

### Standard Features:

- Plastic side panels with 47 insertion opportunities per column
- Key locking roll up tambour door
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand
- Twelve year warranty against manufactures defect

### Assembled cart dimensions:

- 80.5" H x 59.6" W x 27.9" D  
(204.5 cm H x 151.4 cm W x 70.8 cm D)
- Cart net weight: 400 pounds (181 kg)
- Approximate shipping weight: 590 pounds (267.2 kg)