



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Compact Isolation Station
 ISO6548**



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Compact footprint
- Silver ion paint coating with antimicrobial properties
- Unique open well design for visual storage of isolation supplies
- Single folding side shelf with quick-release latch
- Well lid can be latched open
- Built-in glove box holder, holds two glove boxes
- Area for hand sanitizer/soap dispenser
- Adjustable coat rack
- Hard baked paint finish
- Patent pending on design
- Ball-bearing, full extension drawer slides
- 5" Ball bearing premium rubber casters, full swivel, 1 with toe-locking brake and 1 directional
- 15 1/2" high storage compartment with latching door and adjustable shelf
- Cart is delivered fully assembled

Thirteen available color choices:

Standard Color Choice: Cream Yellow
 Call for other choices

Assembled Cart Dimensions

Open: 58" H x 35.75" W x 20" D
 (147.3 cm x 91 cm x 51 cm)

With coat and hat rack extended: 67.4"H

Closed: 46.7" H x 25" W x 20." D
 (118.6 cm x 63.5 cm x 51 cm)

Cart net weight: 98 pounds (44.6 kg)

Approx. shipping weight: 132 pounds (60 kg)

Open Well Dimensions:

- 13"H x 21.375"W x 16" D
 (33 cm H x 54.29 cm W x 40.64 cm D)

Bottom Cabinet Internal Dimensions:

- 16.375"H x 21.375"W x 16" D
 (41.6 cm H x 54.29 cm W x 40.64 cm D)