

MR-Conditional Line/ Model: MR7B-EMG

Aluminum Seven Drawer Anesthesia Cart, Breakaway Lock, Specialty Package



Over 20 Drawer Colors Available:

Antimicrobial Finish Available: Lemon Drop Yellow Only

Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Emergency Specialty Package Includes MRI-Compatible Accessories:

- Stainless steel IV pole (MR-IV2)
- Plastic cardiac board and painted aluminum mounting brackets (MR-CARDBRD)
- Painted aluminum oxygen tank holder – D/E size cylinders (MR-ALO2)
- Aluminum utility hooks (4 ea.) (MR-UHOOKS4)
- Plastic full drawer trays with adjustable dividers (2 ea.) (68530-P1)

Construction:

- Non-magnetic or "weekly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- Painted aluminum breakaway locking system
- Breakaway lock seals, numbered, yellow, pkg-100
- Protective corner bumpers
- Painted aluminum pull-out shelf
- Antimicrobial plastic top with integrated push handles, removable
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (4) - 3.25" (8.3 cm)
- (1) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Internal Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Overall Cart Dimensions:

- 45.75" H x 30.5" W x 23.75" D (116.2 cm H x 77.5 cm W x 60.3 cm D)
- Writing surface: 43.5" (110.5 cm)
- Cart net weight: 127 pounds (57.6 kg)
- Approximate shipping weight: 163 pounds (74.1 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7 cm H x 94 cm W x 66 cm D)