

## PowerSpin™ FX Centrifuge

### Product Images



### Short Description

1. PowerSpin™ FX CENTRIFUGE
2. FIXED SPEED 3400 RPM
3. 6 PLACE FIXED ANGLE ROTOR
4. 30 MIN. MECHANICAL TIMER
5. 6X10ML OR 3X15ML CAPACITY
6. POWER INPUT 110V

## Description

---

UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8- 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.

## Additional Information

Accessories	Header	Accessories
	C800-00	8 Place Rotor
	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster

Specification	Number of Tubes	6
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Fixed
	Rotor	Fixed Angle
	Maximum RPM	3400
Maximum RCF(g)	1320	

Size and Weight	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	13 lbs
	Shipping Weight	15 lbs