

Standard Line 294 Capacity Unit Dose Medication Box Cart

#SL294BOX

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 294 medication box capacity, single width narcotics drawer
- BEST® brand lock on cabinet with four keys
- Illinois lock on narcotics drawer with four keys
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Pull-out side shelf
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Dividers are available to accommodate different medication packaging

Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 6 ea – 3.25" med box drawers with internal dividers, capacity 42 boxes per drawer
- 1 ea – 9.75" multipurpose storage drawer with adjustable dividers

Right Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D
(116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 pounds (95.5 kg)
- Approx. shipping weight: 255 pounds (116 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D
(144.8 cm H x 119.4 cm W x 68.6 cm D)

Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts