



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## ETC Line, Four Drawer, Treatment Cart #ETC-4

### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish



#### Drawer Configuration:

- (1) – 2.75" (6.98 cm)
- (3) – 5.75" (14.6 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Assembled Cart Dimensions

- 27.25" H x 23.38" W x 15.13" D  
(69.2 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. Shipping Weight: 63 pounds (28.7 kg)

\*\* PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS