



M-series M3DS1830K06+MD18-ANS

M-SERIES TALL ANESTHESIA CART Mini Six Drawer Cart w/MD18-ANS Accessory Package

Tall Height, Mini Width, Six Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-ANS Accessories Include:



3-Tier Back Rail (MD18-RBR3)



(3) 12" Side Rails (SIDERAIL12)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



Engineered & Manufactured in the U.S.



Color shown with a Cream body and Navy drawers. More than 20 additional colors are available.

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Overall Cart Dimensions:

- Work Surface Height: 39.625" (100.6 cm)
- Cart net weight: 95 lbs (43.1 kg)
- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 40.75" H x 18" W x 18" D (103.5 cm H x 45.7 cm W x 45.7 cm D)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)