

# S19LWHPLUS2



20.13" x 17.75" x 17.63" (H x W x D)

Compact refrigerator-freezer with NIST calibrated thermometer, hospital grade cord, and front-mounted lock for general purpose use; replaces S19LWHPLUS

## Highlights:

Designed for general purpose use in medical and scientific facilities

Traceable NIST calibrated thermometer provides external readout of current and high/low temperature

Includes hospital grade cord with a 'green dot' plug

## Product Features:

Compact dimensions	Sized for easy use on top of counters
NIST calibrated temperature display	Thermometer provides continuous display of the current and high/low temperature to the nearest tenth of a degree, with a readout switchable from Celsius to Fahrenheit at the touch of a button
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety
Factory installed lock	Keyed lock ensures a secure interior
Manual defrost	Static cooling system for improved temperature stability
Door storage	Store items right on the door for easy access
Adjustable thermostat	Dial thermostat conveniently located on the door to the freezer compartment for easy temperature management
Removable shelf	Wire shelf included for storage convenience

## S19LWHPLUS2 Specifications:

Overview	
Height of Cabinet	19.75" (50 cm)
Height to Hinge Cap	20.13" (51 cm)
Width	17.75" (45 cm)
Depth	17.63" (45 cm)
Depth with door at 90°	33.88" (86 cm)
Capacity	1.7 cu.ft. (48 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.7
Voltage/Frequency	115 V AC/60 Hz
Weight	31.0 lbs. (14 kg)
Shipping Weight	66.0 lbs. (30 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101057873
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Refrigerator - Interior Height	16.13" (41 cm)
Refrigerator - Interior Width	14.88" (38 cm)
Refrigerator - Interior Depth	14.25" (36 cm)
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	1
Refrigerator - Full Door Shelves	1
Refrigerator - Half Door Shelves	1
Freezer - Interior Height	4.38" (11 cm)
Freezer - Interior Width	8.75" (22 cm)
Thermostat Type	Dial
Refrigerant Type	R600a
Refrigerant Amount	0.42 oz.
High Side PSI	120.0

Low Side PSI	70.0
Compressor Step Height	7.38" (19 cm)
Compressor Step Width	14.88" (38 cm)
Compressor Step Depth	6.0" (15 cm)

