

FF28LWHVAC



24.88" x 18.5" x 17.75" (H x W x D)

Compact medical all-refrigerator for temperature stable medical storage, with digital thermostat, lock, self-closing door, temperature alarm, and internal fan

Highlights:

Medical grade design conforms with CDC guidelines for vaccine storage

Specialized interior design protects stored contents when the door is open

Automatic defrost with an internal fan for stabilized cooling designed to stay within 2°C of set temperature

Product Features:

Improved temperature stability	Interior is designed to operate within 2°C of set temperature performance
Interior drawer covers	Two shelves are protected by metal door flaps that help to protect stored contents from being exposed to room temperature when the door is opened
Digital thermostat	Electronic controls ensure precise and intuitive temperature management
High/low temperature alarm	Audible alarm sounds if the interior temperature rises or falls out of set range
NIST calibrated temperature display	Alarm includes NIST calibrated display of the current interior and high/low temperature
Self-closing door	Door closes automatically when opened at angles less than 90° to prevent the interior from being exposed to room temperature
3/8" access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Hospital grade cord and plug	Secure three-pronged 'green dot' cord for added unit safety in high traffic settings
Automatic defrost	Reduced maintenance with auto defrost system

FF28LWHVAC Specifications:

Overview	
Height of Cabinet	24.5" (62 cm)
Height to Hinge Cap	24.88" (63 cm)
Width	18.5" (47 cm)
Depth	17.75" (45 cm)
Depth with door at 90°	34.5" (88 cm)
Capacity	0.8 cu.ft. (23 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.7
Voltage/Frequency	115 V AC/60 Hz
Weight	60.0 lbs. (27 kg)
Shipping Weight	63.0 lbs. (29 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101007175
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.7oz.
High Side PSI	120.0
Low Side PSI	70.0
Level Legs Qty	2
Temperature Range	2 to 8°C
Dimensions	
Interior Height	11.38" (29 cm)
Interior Width	14.0" (36 cm)
Interior Depth	10.5" (27 cm)

