



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Value Line Med Bin Cart with 36 – 3.5 inch Bins, Basic Electronic Pushbutton Lock with Key Lock Override #VLT36EBIN3

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.

### Specifications



Color Shown is Sand



#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Tall cabinet Basic E-lock system, key lock override with 4 keys
- Factory set user and supervisor code
- Lock can be programmed in field
- Internal locking narcotics box in bottom 9" drawer
- Pull out side shelf
- Includes 36 Harloff 3.5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

#### Drawer Configuration:

- 1 – ea 3" multipurpose storage drawers
- 6 – ea 3" bin drawers including 6 - 3.5" bins each
- 1 – ea 9" drawer with 1 divider and locking narcotics box

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

#### Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D  
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (94.5 kg)

#### Thirteen available color choices:

(may be viewed at [www.harloff.com](http://www.harloff.com))

- Beige, Burgundy, Light Blue, Dark Blue, Light Gray, Mauve, Red, Teal, Yellow, Hammertone Blue, Hammertone Gray, Hammertone Green and Sand