



M-series M3DS2430B06+MD24-EMG

M-SERIES TALL EMERGENCY CART with MD24-EMG Emergency Package

Tall, Medium Width, Six Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD24-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board and brackets (MD-CARDBRD)

Color shown is Red. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Latex free
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064

medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300

Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Assembled Cart Dimensions:

- 40.75" H x 30" W x 22" D (103.5 cm H x 76.2 cm W x 55.9 cm D) cart only, 68.8" H x 32.7" W x 22" D (174.7 cm H x 83.1 cm W x 55.9 cm D) w/accessories
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts