

PHC61G



39.5" x 23.63" x 23.5" (H x W x D)

Single chamber warming cabinet with glass door, stainless steel interior, digital thermostat, and lock

Highlights:

Stainless steel construction inside and out for durable use and easy cleaning

Digital controls conveniently located for easy adjustment and temperature viewing

Full length glass door keeps complete interior in view for easy inventory check

Product Features:

Digital thermostat	Electronic controls are conveniently located for easy temperature adjustment
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature
120 to 150°F Temperature Range	Designed and set for hospital blankets
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior
Professional handle	Stainless steel handle completes the professional look
Cantilevered shelving system	Easily adjustable to accommodate virtually any size
Stainless steel interior	Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning
Factory installed lock	Keyed lock for a secure interior
Reversible door	User-reversible door swing for added flexibility

PHC61G Specifications:

Overview

Height of Cabinet	39.5" (100 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Depth with Handle	25.38" (64 cm)
Depth with door at 90°	45.5" (116 cm)
Capacity	5.0 cu.ft. (142 L)
Door	Glass
Cabinet	Black
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Shipping Weight	145.0 lbs. (66 kg)
Estimated Time to Ship	Next day shipping
Parts & Labor Warranty	1 Year

Warming Cabinet

Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Shelf Type	Wire
Shelf Qty	2
Thermostat Type	Digital
Interior Height	26.0" (66 cm)
Interior Width	20.0" (51 cm)
Interior Depth	19.5" (50 cm)
Dolly	Optional
Temperature Range	120 to 150°F

