

AGP34RFLCALADA



32.75" x 19" x 19.75" (H x W x D)

General purpose ADA height 2-door refrigerator-freezer with digital refrigerator control, internal fan, combination lock, and NIST calibrated alarm/thermometer

Highlights:

Combination refrigerator-freezer for general purpose storage in facilities

Digital controls allow you easily manage the refrigerator temperature between 2 and 8°C

NIST calibrated thermometer included to help monitor the interior temperature

Product Features:

Refrigerator-freezer combination	2 to 8°C refrigerator section for basic storage applications, with an upper freezer compartment ideal for ice packs
Combination lock	Factory-installed front lock secures both compartments with a user-set code. Two keys are included to override the combination as needed
Digital thermostat (refrigerator)	Easy operation with a factory-installed digital thermostat that allows you to manage the refrigerator temperature between 2 and 8°F (NOTE: the freezer temperature will self-adjust based on the refrigerator setting)
Alarm/thermometer included	Temperature monitoring device (TMD) displays the current and min/max temperature of the refrigerator interior, with an audible alarm to alert if the unit goes out of range
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature
NIST calibrated external temperature readout	Thermometer is NIST calibrated in our ISO/IEC 17025:2017 certified laboratory and ships with a 2-year certificate of calibration for your records
Fan-forced cooling	Refrigerator section includes a factory-installed fan to help ensure a more evenly cooled interior
Automatic defrost refrigerator	Lower refrigerated section does not require manual defrosting

AGP34RFLCALADA Specifications:

Overview	
Height of Cabinet	32.63" (83 cm)
Height to Hinge Cap	32.75" (83 cm)
Width	19.0" (48 cm)
Depth	19.75" (50 cm)
Depth with door at 90°	36.75" (93 cm)
Capacity	3.2 cu.ft. (91 L)
Defrost Type	Cycle
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	55.0 lbs. (25 kg)
Shipping Weight	77.0 lbs. (35 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086033
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	18.5" (47 cm)
Refrigerator - Interior Width	15.75" (40 cm)
Refrigerator - Interior Depth	15.0" (38 cm)
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	2
Refrigerator - Full Door Shelves	2
Refrigerator - Temperature Range	2 to 8°C
Freezer - Interior Height	7.5" (19 cm)
Freezer - Interior Width	14.0" (36 cm)
Freezer - Interior Depth	14.0" (36 cm)
Freezer - Full Door Shelves	1

Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	1.16 oz.
High Side PSI	250.0
Low Side PSI	88.0
Level Legs	2
Compressor Step Height	9.0" (23 cm)
Compressor Step Width	15.75" (40 cm)
Compressor Step Depth	5.0" (13 cm)

