

# FS30LVAC



24.63" x 18.5" x 17.75" (H x W x D)

Countertop vaccine storage all-freezer with digital thermostat, temperature alarm, hospital grade cord, self-closing door, and front lock

## Highlights:

Medical grade design conforms with CDC guidelines for vaccine storage

Compact size is ideal for countertop use

Digital thermostat ensures more precision in temperature control

## Product Features:

|                                           |                                                                                                                          |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Compact dimensions                        | Sized for easy use on top of counters                                                                                    |
| Digital thermostat                        | Electronic controls ensure easy and intuitive temperature management, with a readout in Celsius or Fahrenheit            |
| High/low temperature alarm                | Audible alarm sounds if the interior temperature rises or falls out of set range                                         |
| NIST calibrated temperature display       | Alarm includes NIST calibrated display of the current interior and high/low temperature to the nearest tenth of a degree |
| Self-closing door                         | Door closes automatically when opened at angles less than 90°                                                            |
| 3/8" access port                          | Unit comes with a factory installed probe hole for users to connect their own monitoring equipment                       |
| Hospital grade cord with 'green dot' plug | Secure three-pronged cord for added unit safety                                                                          |
| Manual defrost                            | Static cooling system for improved temperature stability                                                                 |
| Factory installed lock                    | Keyed lock for a secure interior                                                                                         |

## FS30LVAC Specifications:

| Overview                   |                    |
|----------------------------|--------------------|
| Height of Cabinet          | 26.25" (67 cm)     |
| Height to Hinge Cap        | 24.63" (63 cm)     |
| Width                      | 18.5" (47 cm)      |
| Depth                      | 17.75" (45 cm)     |
| Depth with door at 90°     | 34.5" (88 cm)      |
| Capacity                   | 1.8 cu.ft. (51 L)  |
| Defrost Type               | Manual             |
| Door                       | White              |
| Cabinet                    | White              |
| US Electrical Safety       | ETL                |
| Canadian Electrical Safety | ETL-C              |
| Amps                       | 0.8                |
| Voltage/Frequency          | 115 V AC/60 Hz     |
| Weight                     | 50.0 lbs. (23 kg)  |
| Shipping Weight            | 72.0 lbs. (33 kg)  |
| Parts & Labor Warranty     | 2 Years            |
| Compressor Warranty        | 5 Years            |
| UPC                        | 761101032917       |
| Freezer                    |                    |
| Door Swing                 | RHD                |
| Reversible                 | Factory Reversible |
| Interior Height            | 18.63" (47 cm)     |
| Interior Width             | 13.38" (34 cm)     |
| Interior Depth             | 13.0" (33 cm)      |
| Shelf Type                 | Wire               |
| Shelf Qty                  | 2                  |
| Thermostat Type            | Digital            |
| Refrigerant Type           | R600a              |
| Refrigerant Amount         | 0.8 oz.            |
| Level Legs                 | 2                  |
| Compressor Step Height     | 8.0" (20 cm)       |
| Compressor Step Width      | 13.38" (34 cm)     |
| Compressor Step Depth      | 5.75" (15 cm)      |

