

FF590SSHMED



34.13" x 23.63" x 24.25" (H x W x D)

24" wide frost-free built-in undercounter all-refrigerator with temperature alarm, thermometer, lock, internal fan, green dot cord, and stainless steel door

Highlights:

Flexible design allows built-in or freestanding use in 24" wide spaces

Fully frost-free interior for low maintenance operation

Temperature alarm warns when unit is out of range

Product Features:

Frost-free operation	Eliminates the need to manually defrost the 5.7 cu.ft. interior
Built-in capable	Front-breathing design allows fully recessed installation under counters
Factory installed lock	Keyed lock offers added security
Alarm/thermometer	NIST traceable temperature alarm sounds if the interior goes out of high/low temperature range. Alarm is externally mounted on the door and includes a current and high/low temperature readout
Hospital grade cord	Secure three-pronged cord for added unit safety
Stainless steel door	Professional finish offers lasting durability
Fully finished textured black cabinet	Allows the unit to be used freestanding
Stainless steel handle	Thin (14mm diameter) handle is mounted horizontally for a sleek look
Fan-forced cooling	Interior fan ensures an even temperature throughout the interior

FF590SSHMED Specifications:

Overview	
Height of Cabinet	34.13" (87 cm)
Height to Hinge Cap	34.38" (87 cm)
Width	23.63" (60 cm)
Depth	24.25" (62 cm)
Depth with Handle	25.75" (65 cm)
Depth with door at 90°	47.0" (119 cm)
Capacity	5.7 cu.ft. (161 L)
Defrost Type	Frost-Free
Door	Stainless Steel
Cabinet	Black
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Amps	1.3
Voltage/Frequency	115 V AC/60 Hz
Weight	110.0 lbs. (50 kg)
Shipping Weight	120.0 lbs. (54 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Shelf Type	Glass
Shelf Qty	3
Adjustable Shelves	Yes
Thermostat Type	Dial
Fan Type	Interior and Exterior
Refrigerant Type	R134a
Refrigerant Amount	3.5oz.
High Side PSI	235.0
Low Side PSI	140.0
Interior Light	Yes
Dimensions	
Interior Height	27.63" (70 cm)
Interior Width	21.0" (53 cm)
Interior Depth	17.25" (44 cm)

