



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Three Drawer Mini24 Infection Control Cart (Short Cabinet) with Bumper and 5" Wheels, with Key Lock #4243K

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Yellow.
Additional Colors Available.

Shipping Dimensions:

- 49" H x 33" W x 26" D
(124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 142 pounds (64.4 kg)

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) of drawer space
- Key lock with two keys
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (2) - 9-inch (22.9 cm)

Assembled Cart Dimensions:

- 36.5" H x 29.25" W x 22" D
(92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 103 pounds (46.8 kg)