

FFAR12W



58.38" x 23.63" x 22.75" (H x W x D)

10.1 cu.ft. commercial all-refrigerator with stainless steel interior, white exterior, digital thermostat, lock, and automatic defrost operation

Highlights:

Perfectly sized with a 24" footprint and over 10 cu.ft. of interior storage

Stainless steel interior ensures improved sanitation

Cantilevered shelving allows endless storage configurations to meet your needs

Product Features:

Commercially approved	ETL-S listed to ANSI-NSF Standard 7
Stainless steel interior	Full construction on interior walls and floor ensures improved sanitation
Digital thermostat	Electronic controls for more precise temperature management, with temperature readout and controls externally located for added convenience
Cantilevered shelves	Adjustable chrome shelves can be arranged in virtually any configuration for flexible storage
Factory installed lock	Keyed lock offers convenient security
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety
Flat door liner	Interior design conforms to CDC guidelines to ensure stable conditions for sensitive materials
Automatic defrost	Reduced maintenance with auto defrost system
Internal fan	Fan-forced cooling ensures even temperature distribution from top to bottom

FFAR12W Specifications:

Overview	
Height of Cabinet	58.38" (148 cm)
Width	23.63" (60 cm)
Depth	22.75" (58 cm)
Depth with Handle	24.63" (63 cm)
Depth with door at 90°	45.0" (114 cm)
Capacity	10.1 cu.ft. (286 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Other Energy Info	120 Volts / 60 Hz.
Amps	1.7
Weight	145.0 lbs. (66 kg)
Shipping Weight	150.0 lbs. (68 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101049069
Refrigerator Features	
Door Swing	RHD
Reversible	Yes
Shelf Type	Chrome
Shelf Qty	4
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	1.66oz.
High Side PSI	270.0
Low Side PSI	105.0
Interior Light	Yes
Temperature Range	2° to 6°C
Dimensions	
Interior Height	50.25" (128 cm)
Interior Width	19.63" (50 cm)
Interior Depth	16.75" (43 cm)

Compressor Step Height	4.5" (11 cm)
Compressor Step Width	19.63" (50 cm)
Compressor Step Depth	4.5" (11 cm)

