



M-SERIES MEDIUM WIDTH CART

w/Three Drawers & Electronic Pushbutton Lock

Features:

- · Sturdy, uni-body steel construction
- · Basic electronic pushbutton lock
- Key lock override
- · Factory set user and supervisor code
- · Cabinet may be programmed in the field
- · Durable scratch resistant heat sealed powder coat finish
- · Full extension drawers slides for easy access
- · Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- · 24" of vertical drawer space
- · Cart is delivered assembled
- Industry only fully curved protective bumper
- · Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- · Latex free

Colors:

Shown in Yellow, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Short Height, Medium Width, Three Drawer Cart with **Electronic Pushbutton Lock**

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- One 6" drawer (15.2 cm) and two 9" drawers (22.9 cm)

Assembled Cart Dimensions:

- 37.25" H x 30" W x 22" D (94.6 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
 - Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on

shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at

https://www.harloff.com/product-category/m-series-medical-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300