



# M-series

## MDS2424ELP05

### M-SERIES MEDIUM WIDTH CART w/Five Drawers, Electronic Keypad Lock & Prox

#### Features:

- Sturdy, uni-body steel construction
- Durable scratch resistant heat sealed powder coat finish
- Electronic locking cart with keypad or prox reader access and Key Lock override
- Easily add new users with existing HID compatible badges (Prox and iClass)
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user and user groups
- Battery operated with six D-cell Batteries; life 1 year minimum at 25-50 uses/day
- Lock has an LED low battery indicator
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Engineered &  
Manufactured  
in the U.S.



Color shown in  
Forest Green.  
More than 20  
additional colors  
are available.



#### Short Height, Medium Width, Five Drawer Cart with an Electronic Keypad Lock & Proximity Reader

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)

#### Assembled Cart Dimensions:

- 37.25" H x 30" W x 22" D (94.6 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at  
<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300