



M-series

MDS2430K06

M-SERIES MEDIUM WIDTH CART with Six Drawers and a Key Lock

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
 - 1 steering caster and 2 braking casters
- Latex free



Engineered &
Manufactured
in the U.S.

Color shown is Light Gray with
Hammertone Gray drawers. More than
20 additional colors are available.



Tall, Medium Width, Six Drawer Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.9 cm)

Assembled Cart Dimensions:

- 43.75" H x 30" W x 22" D (111.1 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 119 pounds (54 kg)
- Approx shipping weight: 158 pounds (71.7 kg)
- Approx shipping dims: 56" H x 33" W x 26" D (142.2 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300