



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Value Line Med Bin Cart with 12 – 5 inch Bins, Key Locking
 #VLT12BIN5** **Specifications**



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system
- Internal locking narcotics box in bottom 9” drawer
- Pull out side shelf
- Includes 12 Harloff 5” medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five year warranty
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green

Drawer Configuration:

- 1 - ea 3” multipurpose storage drawer
- 3 - ea 3” bin drawers including 4 - 5” bins each
- 1 - ea 9” drawer with 2 adjustable dividers
- 1 - ea 9” drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 145 pounds (65.9 kg)
- Approx. shipping weight: 180 pounds (81.8 kg)

Value Line Drawers Have the Ability to Easily Convert in the Field:

- If your medication packaging changes, the drawers in your Value Line Cart can change
- Replace 3 small drawers with 1 large drawer



Unit Dose/Multi Dose Drawer



Punch Card Drawer