

ACR1151



58.38" x 23.63" x 22.88" (H x W x D)

Mid-sized pharmaceutical all-refrigerator with stainless steel construction inside and out, digital controls, and self-closing glass door

Highlights:

Stainless steel construction inside and out for easier cleaning and improved sanitation

Self-closing door and open door alarm ensure items are not accidentally overexposed to warm temperatures

Full display all the way down to the floor

Product Features:

Stainless steel interior and exterior finish	Cabinet and interior walls are constructed from stainless steel
Frost-free	No-frost operation ensures minimum user maintenance
Self-closing door	Door will automatically close when left ajar, preventing items from accidentally being exposed to ambient temperatures
Glass door	Tempered glass door with stainless steel trim provides heat-safe view of interior contents
Recessed LED lighting	Energy efficient lighting offers full view of interior
Factory installed lock	Keyed lock for a secure interior
Adjustable cantilevered shelves	Adjustable shelving accommodates virtually any size with convenient cantilevered system
Fan-cooled interior	Internal fan ensures the interior is evenly cooled from top to bottom
Operates between 2° and 8°C	Designed to maintain temperatures as recommended by the CDC

ACR1151 Specifications:

Overview	
Height of Cabinet	58.38" (148 cm)
Width	23.63" (60 cm)
Depth	22.88" (58 cm)
Depth with Handle	24.0" (61 cm)
Depth with door at 90°	45.0" (114 cm)
Capacity	11.0 cu.ft. (311 L)
Defrost Type	Frost-Free
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.7
Voltage/Frequency	115 V AC/60 Hz
Weight	141.0 lbs. (64 kg)
Shipping Weight	146.0 lbs. (66 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101050348
Refrigerator Features	
Door Swing	RHD
Reversible	No
Shelf Type	Chrome
Shelf Qty	6
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior and Exterior
Refrigerant Type	R600a
Refrigerant Amount	1.66oz.
High Side PSI	270.0
Low Side PSI	105.0
Level Legs Qty	4
Interior Light	Yes
Temperature Range	2° to 8°C
Dimensions	
Interior Height	50.25" (128 cm)
Interior Width	19.38" (49 cm)
Interior Depth	17.75" (45 cm)

Compressor Step Height	4.75" (12 cm)
Compressor Step Width	19.38" (49 cm)
Compressor Step Depth	5.75" (15 cm)

