



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Eight Bronchoscope/ENT Scope Storage Cabinet, with Tambour Locking Door, Standard Package #SCBR-8

The SCBR-8 is a state of the art stationary cabinet that stores and dries up to 8 Broncho/ENT scopes.

Specifications



**Color Show is Sand
Custom Colors Available**

Construction:

- Sturdy, 18 gauge steel

Features:

- Holds up to 8 scopes in a secure stationary location
- Back wall with padded vinyl prevents damage to scope optics
- Lockable tambour door helps protect and secure scopes
- Engineered to protect scopes and provide easy access
- Removable drip tray

Drying Package Available: #SCBR-8DP Includes

- Drying fans 100-240VAC Universal Input (220 with International plug adapter # SCINTLPLUGADPTR)
- HEPA filter captures up to 99.97% of particles, 0.3 microns, exceeding AORN scope storage recommendations

Assembled Cart Dimensions:

- 56.125"H x 36.125"W x 18.75"D
(142.6cm H x 91.8 cm W x 47.6 cm D)
- Cart net weight 175 pounds (79.4 kg)