



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

V-Series, Six Drawer Crash Cart with Breakaway Lock #V30-5B

This high capacity crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price.

Specifications

Construction:

- Steel construction

Features:

- Break-away lock bar
- 100 break-away seals included
- 30" of vertical drawer space
- Aluminum side rail / push handle
- Corner bumpers, protecting the cart
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)!

Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (2) – 9" (30.5 cm)

Integrated Pull Handle on all drawers

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D)
- Writing surface: 40" (101.6cm)
- Cart weight: 117 pounds (53.2 kg)
- Approx. shipping weight: 158 pounds (71.8 kg)
- Approx. shipping dims: 57" h x 33" W x 26" D (144.8cm H x 83.82cm W x 66cm D)



Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black