


FF28LWHPRO



25" x 18.5" x 17.75" (H x W x D)

Compact auto defrost all-refrigerator with digital thermostat, internal fan, front lock, and probe hole for user-installed monitoring equipment

Highlights:

Advanced precision with digital temperature controls

Factory installed access port allows users to install their own monitoring equipment

Compact size is ideal for use on top of counters

Product Features:

Digital thermostat	Digital microprocessor for easy and precise temperature management, with all controls externally located for added convenience
Factory installed lock	Keyed lock for added security, with two keys included
Automatic defrost	Reduced maintenance with auto defrost system
Internal fan	Fan-forced cooling ensures the unit is evenly cooled from top to bottom
Flat door liner	Interior design conforms to CDC guidelines to ensure stable conditions for sensitive materials
3/8" access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Adjustable shelving	Includes one full-width adjustable wire shelf and one removable half shelf that can be positioned above the compressor step
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety in high traffic areas
Foamed-in-place cabinet	Quality design ensures efficiently insulated interior

FF28LWHPRO Specifications:

Overview	
Height of Cabinet	25.0" (64 cm)
Width	18.5" (47 cm)
Depth	17.75" (45 cm)
Depth with door at 90°	34.5" (88 cm)
Capacity	2.4 cu.ft. (68 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.7
Voltage/Frequency	115 V AC/60 Hz
Weight	45.0 lbs. (20 kg)
Shipping Weight	48.0 lbs. (22 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101057231
Refrigerator Features	
Door Swing	RHD
Reversible	Yes
Shelf Type	Wire
Shelf Qty	2
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.7oz.
High Side PSI	120.0
Low Side PSI	70.0
Level Legs Qty	2
Temperature Range	2 to 8°C / 36 to 46°F
Dimensions	
Interior Height	21.0" (53 cm)
Interior Width	15.5" (39 cm)
Interior Depth	13.0" (33 cm)
Compressor Step Height	7.75" (20 cm)

Compressor Step Width

Compressor Step Depth

