

FFAR10PRO



54.88" x 23.63" x 23" (H x W x D)

24" wide 10 cu.ft. auto defrost commercial all-refrigerator with lock, digital thermostat, internal fan, and access port for user-provided monitoring equipment

Highlights:

Factory installed access port allows users to install their own monitoring equipment

Advanced precision with digital temperature controls

Ideal 24" footprint offers over 10 cu.ft. of storage in a slim fit

Product Features:

Large capacity	10.1 cu.ft. capacity in a 24" footprint
Automatic defrost	Reduced maintenance with auto defrost system
Digital thermostat	Electronic controls ensure easy and intuitive temperature management
Access port	Pre-installed 3/8" probe hole included in the rear allows users to connect their own monitoring equipment
Factory installed lock	Top-mounted for convenient security (NOTE: lock can alternatively be installed on the unit's side. Specify upon ordering)
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety
Internal fan	Fan-forced cooling ensures the unit is evenly cooled from top to bottom
Adjustable shelves	Rearrange refrigerator space to accommodate all sizes or remove shelves for simple cleanup
Flat door liner	Easy to clean, with more depth for storage

FFAR10PRO Specifications:

Overview	
Height of Cabinet	56.0" (142 cm)
Height to Hinge Cap	54.88" (139 cm)
Width	23.63" (60 cm)
Depth	23.0" (58 cm)
Depth with Handle	24.75" (63 cm)
Depth with door at 90°	45.38" (115 cm)
Capacity	10.1 cu.ft. (286 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Other Energy Info	120 Volts / 60 Hz.
Amps	1.5
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Shipping Weight	145.0 lbs. (66 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101057194
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Shelf Type	Wire
Shelf Qty	4
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R134a
Refrigerant Amount	3.0oz.
High Side PSI	290.0
Low Side PSI	70.0
Level Legs Qty	2
Interior Light	Yes
Temperature Range	2 to 8°C / 36 to 46°F
Dimensions	
Interior Height	48.0" (122 cm)

Interior Width	21.0" (53 cm)
Interior Depth	17.75" (45 cm)
Compressor Step Height	6.0" (15 cm)
Compressor Step Width	21.0" (53 cm)
Compressor Step Depth	5.5" (14 cm)

