

# ACR1151-FS24LSTACKPRO



78.5" x 24" x 23" (H x W x D)

Compact manual defrost FS24LPRO all-freezer stacked with 11 cu.ft. pharmaceutical refrigerator ACR1151PRO, both with factory-installed probe holes

### Highlights:

Pharmaceutical all-refrigerator features complete stainless steel construction inside and out for lasting durability

Separate refrigerator and freezer cooling in one slim-fitting footprint

Both models are equipped with factory-installed probe holes for user-provided monitoring equipment

### Product Features:

FS24LPRO freezer	Compact manual defrost all-freezer with flat door liner and wire shelf
Access ports	Refrigerator and freezer include a 3/8" factory-installed probe hole to allow users to connect their own monitoring equipment
ACR1151PRO all-refrigerator	Pharmaceutical refrigerator with self-closing glass door, LED lighting, and stainless steel construction inside and out
Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
Factory installed locks	Both units are equipped with front locks, with two keys included for each unit
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Anti-tip bracket included	StackRack ships with anti-tip bracket for easy aftermarket installation
Foamed-in-place cabinets	Designed to ensure a sufficiently insulated interior
Magnetic door gaskets	Quality door gaskets ensure a tight seal for energy savings and improved performance

## ACR1151-FS24LSTACKPRO Specifications:

Overview	
Height of Cabinet	78.5" (199 cm)
Width	24.0" (61 cm)
Depth	23.0" (58 cm)
Depth with door at 90°	45.0" (114 cm)
Capacity	12.4 cu.ft. (351 L)
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Voltage/Frequency	115 V AC/60 Hz
Shipping Weight	215.0 lbs. (98 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101058719

Refrigerator-Freezer	
Door Swing	RHD
Reversible	No
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	50.25" (128 cm)
Refrigerator - Interior Width	19.38" (49 cm)
Refrigerator - Interior Depth	17.75" (45 cm)
Refrigerator - Capacity	11.0 cu.ft
Refrigerator - Shelf Type	Chrome
Refrigerator - Shelf Qty	6
Refrigerator - Temperature Range	2° to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	1.66 oz.
Freezer - Interior Height	14.38" (37 cm)
Freezer - Interior Width	13.5" (34 cm)
Freezer - Interior Depth	13.25" (34 cm)
Freezer - Capacity	1.4 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	1

Thermostat Type	Digital + Dial
Fan Type	Interior
Freezer - Temperature Range	-20 to -15°C
Freezer - Refrigerant Type	R600a
Freezer - Refrigerant Amount	0.8 oz.

