

FF7LW-FS24LSTACKPRO



54" x 23.63" x 23.5" (H x W x D)

FF7LWPRO auto defrost all-refrigerator with digital controls and compact manual defrost FS24LPRO all-freezer with stacking rack, both with factory-installed probe holes

Highlights:

Separate refrigerator and freezer cooling in one 24" wide unit

Both models are equipped with factory-installed probe holes for user-provided monitoring equipment

Customizable design can accept additional options, including temperature alarms, data loggers, and thermometers

Product Features:

Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
FF7LPRO all-refrigerator	Auto defrost all-refrigerator made in Europe and equipped with a digital thermostat and internal fan for optimum temperature performance
FS24LPRO freezer	Compact manual defrost all-freezer with factory-installed probe hole
Access ports	Refrigerator and freezer include a 3/8" factory-installed probe hole to allow users to connect their own monitoring equipment
Hospital grade cords with 'green dot' plugs	Secure three-pronged cords for added unit safety
Factory installed locks	Both models include front locks for added security, with two keys included for each unit
Flat door liner	No door shelving to offer more storage depth and better temperature stability
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Anti-tip bracket included	StackRack ships with anti-tip bracket for easy aftermarket installation

FF7LW-FS24LSTACKPRO Specifications:

Overview	
Height of Cabinet	54.0" (137 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Capacity	6.9 cu.ft. (195 L)
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Amps	1.4
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years

Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	27.75" (70 cm)
Refrigerator - Interior Width	20.0" (51 cm)
Refrigerator - Interior Depth	17.5" (44 cm)
Refrigerator - Capacity	5.5 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	3
Refrigerator - Temperature Range	2 to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	0.64 oz.
Freezer - Interior Height	14.38" (37 cm)
Freezer - Interior Width	13.5" (34 cm)
Freezer - Interior Depth	13.25" (34 cm)
Freezer - Capacity	1.4 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	1
Thermostat Type	Digital + Dial
Fan Type	Interior

Refrigerant Type	R600a
Refrigerant Amount	0.64 oz.
High Side PSI	180.0
Low Side PSI	80.0
Freezer - Temperature Range	-20° to -15°C
Freezer - Refrigerant Type	R134a
Freezer - Refrigerant Amount	1.5 oz.

