

FF7LW-VT65MLSTACKPRO



67" x 23.63" x 23.5" (H x W x D)

FF7LWPRO auto defrost all-refrigerator with digital controls stacked with -25°C manual defrost VT65MLPRO all-freezer, both with factory-installed probe holes

Highlights:

Both models are equipped with factory-installed probe holes for user-provided monitoring equipment

Customizable design can accept additional options, including temperature alarms, data loggers, and thermometers

Separate refrigerator and freezer cooling in one slim-fitting unit

Product Features:

Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
FF7LPRO all-refrigerator	Auto defrost all-refrigerator made in Europe and equipped with a digital thermostat and internal fan for optimum temperature performance
VT65MLPRO all-freezer	Manual defrost -25°C all-freezer with pull-out basket drawers
Factory installed locks	Both models include front locks for added security, with two keys included for each unit
Access ports	Refrigerator and freezer include a 3/8" factory-installed probe hole to allow users to connect their own monitoring equipment
Hospital grade cords with 'green dot' plugs	Secure three-pronged cords for added unit safety
Flat door liner	No door shelving to offer more storage depth and better temperature stability
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Anti-tip bracket included	StackRack ships with anti-tip bracket for easy aftermarket installation

FF7LW-VT65MLSTACKPRO Specifications:

Overview	
Height of Cabinet	67.0" (170 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Capacity	8.7 cu.ft. (246 L)
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Amps	1.4
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years

Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Interior Drawers	Yes
Drawer Qty	3
Interior Light	Yes
Refrigerator - Interior Height	27.75" (70 cm)
Refrigerator - Interior Width	20.0" (51 cm)
Refrigerator - Interior Depth	17.75" (45 cm)
Refrigerator - Shelf Type	Glass
Refrigerator - Shelf Qty	3
Refrigerator - Full Door Shelves	3
Refrigerator - Temperature Range	2 to 8°C / 36 to 46°F
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	0.64 oz.
Freezer - Interior Height	24.25" (62 cm)
Freezer - Interior Width	17.5" (44 cm)
Freezer - Interior Depth	16.5" (42 cm)
Freezer - Capacity	3.5 cu.ft

Thermostat Type	Digital + Dial
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.64 oz.
High Side PSI	180.0
Low Side PSI	80.0
Freezer - Temperature Range	-25° to -15°C
Freezer - Refrigerant Type	R134a
Freezer - Refrigerant Amount	2.76 oz.

