

FFAR10-FS30LSTACKPRO



80" x 23.63" x 23.5" (H x W x D)

FFAR10PRO 10.1 cu.ft. auto defrost all-refrigerator with digital controls and compact manual defrost FS30LPRO all-freezer with stacking rack, both with factory-installed probe holes

Highlights:

Separate refrigerator and freezer cooling in one slim-fitting footprint

Both models are equipped with factory-installed probe holes for user-provided monitoring equipment

Customizable design can accept additional options, including temperature alarms, data loggers, and thermometers

Product Features:

Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
FFAR10PRO all-refrigerator	10.1 cu.ft. auto defrost all-refrigerator with digital controls, internal fan, flat door liner, and adjustable shelving
FS30LPRO freezer	Compact manual defrost all-freezer with factory-installed probe hole
Factory installed locks	Both models include locks (front-mounted on the freezer, side-mounted on the refrigerator) for added security, with two keys included for each unit
Access ports	Refrigerator and freezer include a 3/8" factory-installed probe hole to allow users to connect their own monitoring equipment
Hospital grade cords with 'green dot' plugs	Secure three-pronged cord for added unit safety (both units)
Flat door liner	No door shelving to offer more storage depth and better temperature stability
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Anti-tip bracket included	StackRack ships with anti-tip bracket for easy aftermarket installation

FFAR10-FS30LSTACKPRO Specifications:

Overview	
Height of Cabinet	80.0" (203 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Depth with door at 90°	45.5" (116 cm)
Capacity	12.1 cu.ft. (343 L)
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Voltage/Frequency	115 V AC/60 Hz
Shipping Weight	198.0 lbs. (90 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101058863

Refrigerator-Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	48.0" (122 cm)
Refrigerator - Interior Width	21.0" (53 cm)
Refrigerator - Interior Depth	17.75" (45 cm)
Refrigerator - Capacity	10.1 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	4
Refrigerator - Temperature Range	2 to 8°C / 36 to 46°F
Refrigerator - Refrigerant Type	R134a
Refrigerator - Refrigerant Amount	3.0 oz.
Freezer - Interior Height	18.63" (47 cm)
Freezer - Interior Width	13.38" (34 cm)
Freezer - Interior Depth	13.0" (33 cm)
Freezer - Capacity	1.8 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	2

Thermostat Type	Digital + Dial
Fan Type	Interior
Freezer - Temperature Range	-20° to -15°C
Freezer - Refrigerant Type	R600a
Freezer - Refrigerant Amount	0.8 oz.

