

AFS49ML



83.75" x 55.25" x 31" (H x W x D)

Performance series Pharma-Lab 49 cu.ft. all-freezer in stainless steel

Highlights:

Improved temperature accuracy with a separate temperature probe for continuous product readout

Additional sensor reads interior air temperature for superior temperature control

Adjustable operating control range from -30 °C to -15°C

Product Features:

| | |
|---|---|
| Purpose-built design | Designed and engineered to store vaccines and other sensitive materials under stable temperature conditions |
| Adjustable temperature | Interior is pre-set to -25°C, with an adjustable control range from -30 to -15°C |
| Temperature alarm | Visual and audible alarms alert if the interior temperature rises or falls out of set range |
| Complete stainless steel construction | Interior and exterior are entirely constructed from stainless steel for lasting durability and easy cleanup |
| Forced air cooling | Top-mounted commercial fan ensures rapid temperature recovery |
| Automatic defrost with adjustable cycle defrost | Adjustable cycle defrost operation ensures the unit avoid temperature spikes and ice buildup |
| Self-closing doors | Closes to save energy but will remain open when the door is opened at a 90° angle for easier loading |
| Open door alarm | Alarm will beep 3 times if the door is left open for 30 seconds, 5 times at 60 seconds, and will sound continuously if the door is open for more than 5 minutes |
| 6 inch casters | Ships with a set of casters (2 locking) for optional user installation |

AFS49ML Specifications:

| Overview | |
|----------------------------|----------------------|
| Height of Cabinet | 83.75" (213 cm) |
| Width | 55.25" (140 cm) |
| Depth | 31.0" (79 cm) |
| Depth with Handle | 32.75" (83 cm) |
| Depth with door at 90° | 56.5" (144 cm) |
| Capacity | 49.0 cu.ft. (1388 L) |
| Defrost Type | Frost-Free |
| Door | Stainless Steel |
| Cabinet | Stainless Steel |
| US Electrical Safety | UL |
| Canadian Electrical Safety | ULC |
| Sanitation | ETL-S |
| Amps | 10.5 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 428.0 lbs. (194 kg) |
| Shipping Weight | 450.0 lbs. (204 kg) |
| Estimated Time to Ship | 1 to 2 days |
| Parts & Labor Warranty | 2 Years |
| Compressor Warranty | 5 Years |
| Freezer | |
| Door Swing | French |
| Interior Height | 51.5" (131 cm) |
| Interior Width | 51.5" (131 cm) |
| Interior Depth | 27.0" (69 cm) |
| Shelf Type | Wire |
| Shelf Qty | 8 |
| Adjustable Shelves | Yes |
| Thermostat Type | Digital |
| Interior Light | Yes |
| Fan Type | Interior |
| Refrigerant Type | R290a |
| Refrigerant Amount | 4.6 oz. |
| Wheels | 4 |
| Temperature Range | -30 to -15°C |

