



# M-series

## M3DS2430B06

### M-SERIES TALL EMERGENCY CART with Six Drawers and a Breakaway Lock

#### Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free

#### Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered &  
Manufactured  
in the U.S.

**Color shown is Red. More than 20 additional colors are available.**



### Tall, Medium Width, Six Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and two 9" drawers (22.7 cm)

#### Assembled Cart Dimensions:

- 40.75" H x 30" W x 22" D (103.5 cm H x 76.2 cm W x 55.9 cm D) cart only,
- 68.8" H x 32.7" W x 22" D (174.7 cm H x 83.1 cm W x 55.9 cm D) w/accessories
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300